

Phospholipid characterization by a TQ-MS data based identification scheme

ASMS 2017 MP-406

Tsuyoshi Nakanishi ¹, Masaki Yamada ¹,
Ningombam Sanjib Meitei ^{2,3}

¹ Shimadzu Corporation, Kyoto, Japan,

² PREMIER Biosoft, CA, USA,

³ PREMIER Biosoft, Indore, India

Phospholipid characterization by a TQ-MS data based identification scheme

Overview

Phospholipids play a wide variety of roles such as the structural units of the cellular membrane, the storage and transport of lipids, the precursor of the lipid mediators and the energy source *in vivo*. To deeply understand the biological function of phospholipids, quantitative analysis for the phospholipid composition ratio in the blood and tissues is required. Phospholipids species have many types

of combination of its characteristic head groups and fatty acid composition. Although identification of phospholipid species typically is performed for the data acquired by a high resolution mass spectrometer, in this study, we report phospholipids database search and quantitative analysis by coupling LC-TQ MS with SimLipid software.

Introduction

Phospholipids constitute a lipid bilayer in a cellular membrane in a living cell and the number of phospholipid species reaches up to several thousands due to the combination of its characteristic head groups and its fatty acid composition. It has been also reported that change of phospholipid levels in blood and tissues correlates with various diseases. Therefore, quantitative analysis for phospholipids has been often performed to evaluate the

change of phospholipid species at a concentration. In this study, the carbon nano tube (CNT) probe, known as a long wavelength fluorescent probe, was administered to mice to evaluate the changes in phospholipid in a liver tissue. To analyze the fluctuation of phospholipids, a precursor ion scan/neutral loss scan was performed by a LC-TQ MS. The change of identified phospholipid levels was compared between control and a probe dosing mouse.

Methods

Carbon nano tube (CNT) probe is known as a fluorescent probe for a long-wavelength. This probe was administered to a mouse at a concentration of 300 mg/mL by a tail vein injection (100 μ L). After 5hr of administering, liver tissues were isolated from a control mouse and a ministered mouse under anesthesia. The isolated tissues were rapidly frozen in liquid nitrogen and crushed to some blocks of an appropriate size. Then these tissue blocks were weighed. Furthermore, after crushing frozen tissue blocks by a bead type crusher, phospholipids were extracted by Bligh & Dyer method. Organic phase was recovered and then evaporated. The sample was dissolved with a solution of

$\text{CHCl}_3/\text{MeOH}$ (1:1). Phospholipid profiling by precursor ion scan (PIS) and neutral loss scan (NLS) with LCMS-8060 were executed for the sample diluted with MeOH. In this measurement, phospholipid was analyzed by PIS at m/z 184 focusing on the characteristic head groups of phosphatidylcholine (PC) and sphingomyeline (SM) or NLS at m/z 141 for ethanolamine of phosphatidylethanolamine (PE). The candidate of phospholipids was estimated for each peak detected by PIS and NLS analysis as a result of database search by SimLipid software from PREMIER Biosoft, CA, USA (www.premierbiosoft.com).

Phospholipid characterization by a TQ-MS data based identification scheme

Table 1 Analytical conditions

HPLC condition	
Column	: Phenomenex Kinetex C8 (150 X 2.1 mm, 2.6 mm)
Mobile phase A	: 20 mM Ammonium formate - water
Mobile phase B	: Acetonitrile/2-Propanol (1:1)
Time program (B%), Curved gradient	: 20% (0 min) - 20% (1.0 min) - 40% (2.0 min) - 92.5% (25 min) - 100% (26.0-35.0 min)
Flow rate	: 0.3 mL/min
Injection volume	: 3 mL
Column oven	: 45°C
MS condition (LCMS-8060)	
Ionization	: ESI(+)/(-)
Nebulizer gas flow rate	: 3 L/min
Heating gas flow rate	: 10 L/min
Drying gas flow rate	: 10 L/min
Interface voltage	: 4 kV(+)/-3 kV(-)
Interface temp.	: 300°C
DL temp.	: 250°C
Block heater temp.	: 400°C




Figure 1 PIS for PC and SM at m/z 184 and NLS for PE of $[M+H-141]^+$

Phospholipid characterization by a TQ-MS data based identification scheme

Results

PIS and NLS measurements to detect phospholipid class of PC, SM and PE

To detect phospholipid species in tissue extracts, we performed PIS at m/z 184 as a product ion for PC / SM and NLS of $[M+H-141]^+$ for PE. Both measurements were performed in positive mode. A typical structural formula for PC, SM and PE are shown in Figure 1. As shown in

Figure 1, PIS and NLS were measured by utilizing a loss of its characteristic head group on MS/MS. Figure 2 shows typical mass chromatograms of PIS and NLS for liver tissue extracts. We can see some detected peaks corresponding to PC, SM and PE.




Figure 2 MS chromatograms of PIS for PC and SM at m/z 184 and NLS for PE of $[M+H-141]^+$

In Figure 2, the change of some phospholipids, especially PC and SM can be also confirmed between control and administration of CNT sample. To evaluate the change of phospholipids at a concentration, we first performed a database search to identify the detected phospholipid species. A database search and quantitative analysis were performed by SimLipid software. Scheme 1 shows a workflow of phospholipid profiling by SimLipid software. TQ-MS data (.lcd file) is imported into the program and a database search is carried out on the basis of a observed

mass value of precursor ion/product ion on PIS/NLS measurement. SimLipid software has a theoretical mass value of total product ion/neutral loss masses corresponding to >5,000 species in the TQ database. In this study, a database search was executed at a range of 0.2 Da as a tolerance. As a result, the candidate phospholipid species were estimated for the detected peaks and the intensity for the detected peaks in a group of control or administration of CNT was integrated for quantitative analysis.

Phospholipid characterization by a TQ-MS data based identification scheme

Next, Figure 4 shows the graph plotted by the peak intensity (Total Abundant) of each phospholipid (PC, SM and PE) integrated among each sample group (a normal and a probe administration, n=3). Phospholipid species for PC, SM and PE, which have been detected in all samples, are shown in this graph. These peak intensity were normalized by a value of the tissue weight. An increase of phospholipids which was considered to be the influence of probe administration, was confirmed as Figure 4 showed. In particular, it was confirmed notably

in sphingophospholipids such as SM(38:1), SM(40:3) and SM(42:3). In addition, an increase of phospholipids such as PC(38:6), PC(40:6), PE(38:6) and PE(40:6) which were considered to contain polyunsaturated fatty acids, was observed as well. On the other hand, some phospholipids, PC(34:1) and PE(34:1) have reduced after probe dosing (5hr). These results suggest that increase and decrease in each class of PLs including the same fatty acid composition are correlated.





Figure 4 Changes of phospholipids in liver tissues from a normal and a fluorescent probe administration

Phospholipid characterization by a TQ-MS data based identification scheme

Conclusions

- We performed phospholipid profiling by coupling TQ-MS measurement with SimLipid database search.
- To evaluate the phospholipid profiling system, we analyzed the change of phospholipids in liver tissues between control and probe administration.

First Edition: June, 2017



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.