

“Like New” Sale on MPC-2200 Multipurpose Sample Compartment Chamber

50% Off List Price

The MPC-2200 enables both reflectance and transmittance measurement of samples having a wide variety of shapes. An integrating sphere is built-in to permit accurate measurement of solid samples. The sample space around the integrating sphere is ample enough to allow measurement of very large samples.

Key features include:

- Measurement wavelength range: 240 to 800 nm
- Maximum sample size:
 - Transmittance: 305 mm dia. × 50 mm thick, or
204 mm dia. × 300 mm thick
 - Reflectance: 305 mm dia. × 50 mm thick
- With independent S/R beam switching, 0°/8° incidence angle reflectance measurement is possible without tilting the sample
- Integrating sphere shift function expands application range
- V stage built-in; sample position can be adjusted vertically and laterally



Please reference the following:

Part number: 206-23830-91
Serial number: A106546001-SITC2-7427

Contains all parts and manuals, fully assembled, no dents, dings, or scratches. Box shows signs of shipping.

Valid in USA only.

Not valid with any other offer. Available on a first-come, first-served basis.

Shimadzu Scientific Instruments, 7102 Riverwood Drive, Columbia, MD 21046
Phone: 800-477-1227; Fax: 410-381-1222; webmaster@shimadzu.com; www.ssi.shimadzu.com